


LifeSpan[®] ePTFE Vascular Grafts

EXCELLENT HANDLING AND SUTURABILITY

 LeMaitre[®]

40 YEARS
1983-2023



LifeSpan® ePTFE Vascular Grafts

EXCELLENT HANDLING AND SUTURABILITY

LifeSpan ePTFE Vascular Grafts for peripheral vascular reconstructions and AV Access have been used in clinical practice for more than 25 years. It features soft ePTFE material with superb suture retention, high burst strength and excellent handling.

- Peripheral and AV Access
- Over 25 years in Clinical Use
- Excellent Suture Retention
- High Burst Strength
- Color Coded Packaging

LIFESPAN GRAFTS FOR AV ACCESS

Offering various sizes, the robust regular wall thickness makes it an ideal implant. Center spiral grafts are recommended for procedures requiring enhanced resistance to kinking and compression. Stepped and Quick Tapered models have been designed to reduce the risk of steal syndrome and high cardiac output.

LIFESPAN GRAFT TECHNICAL DATA

Material

Vascular Graft: ePTFE - expanded PolyTetraFluoroEthylene

External Spiral Support: PTFE - PolyTetraFluoroEthylene

Wall Thickness (6 mm graft - nominal)


Regular/Standard Wall:	0.63 mm
Thin Wall:	0.4 mm

Internodal Distance (nominal): 20 +/- 10 µm

Suture Retention Strength: min. 300 grams

Average (+/- SD) Min. Burst Strength*: 218 +/- 31 psi





* Data on file at LeMaitre Vascular, Inc.



LIFESPAN GRAFTS FOR PERIPHERAL VASCULAR PROCEDURES

Offering regular and thin wall options for reconstructions, the external spiral support can be removed, facilitating an easy creation of the anastomosis to the vessel. Models with full or partial external spiral support are recommended for extra-anatomic surgical reconstructions and those requiring enhanced resistance to kinking and compression.

ORDERING INFORMATION

	Diameter	Length	Spiral Support	Color	Model #	
Regular Wall						
Straight 	6 mm	10 cm		Yellow	R06010	
	6 mm	20 cm		Yellow	R06020	
	7 mm	20 cm		Yellow	R07020	
	8 mm	20 cm		Yellow	R08020	
	6 mm	30 cm		Yellow	R06030	
	6 mm	50 cm		Yellow	R06050	
	7 mm	50 cm		Yellow	R07050	
	8 mm	50 cm		Yellow	R08050	
	6 mm	80 cm		Yellow	R06080	
	7 mm	80 cm		Yellow	R07080	
8 mm	80 cm		Yellow	R08080		
Thin Wall						
Straight 	6 mm	20 cm		Red	T06020	
	6 mm	50 cm		Red	T06050	
	7 mm	50 cm		Red	T07050	
	8 mm	50 cm		Red	T08050	
	6 mm	80 cm		Red	T06080	
	7 mm	80 cm		Red	T07080	
	8 mm	80 cm		Red	T08080	
	10 mm	80 cm		Red	T10080	
	Regular Wall - External Spiral Support					
	Straight 	6 mm	50 cm	50 cm	Yellow	R06050C50
7 mm		50 cm	50 cm	Yellow	R07050C50	
8 mm		50 cm	50 cm	Yellow	R08050C50	
6 mm		80 cm	80 cm	Yellow	R06080C80	
7 mm		80 cm	80 cm	Yellow	R07080C80	
8 mm		80 cm	80 cm	Yellow	R08080C80	
10 mm	80 cm	80 cm	Yellow	R10080C80		
Thin Wall - External Spiral Support						
Straight 	6 mm	50 cm	50 cm	Red	T06050C50	
	7 mm	50 cm	50 cm	Red	T07050C50	
	8 mm	50 cm	50 cm	Red	T08050C50	
	6 mm	80 cm	80 cm	Red	T06080C80	
	7 mm	80 cm	80 cm	Red	T07080C80	
	8 mm	80 cm	80 cm	Red	T08080C80	
	10 mm	80 cm	80 cm	Red	T10080C80	

ORDERING INFORMATION

LifeSpan[®] ePTFE Vascular Grafts

	Diameter	Length	Spiral Support	Color	Model #
Regular Wall - Center Spiral Support					
Straight 	6 mm	50 cm	10 cm		R06050CS
Thin Wall - Center Spiral Support					
Straight 	6 mm	50 cm	30 cm		T06050C30
	8 mm	50 cm	30 cm		T08050C30
	6 mm	80 cm	50 cm		T06080C50
	7 mm	80 cm	50 cm		T07080C50
	8 mm	80 cm	50 cm		T08080C50
Regular Wall					
Stepped 	4-7 mm	50 cm			RS47050
Regular Wall - Center Spiral Support					
Stepped 	4-7 mm	50 cm	10 cm		RS47050CS
Regular Wall					
Quick Tapered 	4-7 mm	40 cm			QT47040
	4-7 mm	50 cm			QT47050
Regular Wall - Center Spiral Support					
Quick Tapered 	4-7 mm	50 cm	10 cm		QT47045CS

These specifications are not intended as a warranty. In the interest of product improvement, these specifications may be changed from time to time without notice. Please consult your sales representative for details.



LeMaitre and LifeSpan are registered trademarks of LeMaitre Vascular, Inc. ©2023 LeMaitre Vascular, Inc. All rights reserved. Printed in the USA. M0275 Rev. J 01/23

LeMaitre Vascular, Inc.
63 Second Avenue
Burlington, MA 01803
USA

T 781 221 2266
F 781 221 2223

LeMaitre Vascular ULC
9135 Keele Street, Suite B6
Vaughan, ON L4K 0J4
Canada

T 905 832 8077
F 905 303 3903

LeMaitre Vascular GmbH
Otto-Volger-Str. 5a/b
65843 Sulzbach/Ts.
Germany

T +49 (0)6196 659230
F +49 (0)6196 527072

LeMaitre Vascular, Pte. Ltd.
138 Robinson Road
#12-02 Oxley Tower
068906
Singapore

T +65 6813 0281

Customer Service
US 800 628 9470
Canada 855 673 2266

www.lemaitre.com