SPHERA DR MRITM SURESCANTM

PACEMAKERS (IPG)

Completely Automatic - Simple to Use

Automatically monitors for implant detection and continuously adapts key device parameters to ensure the therapies are optimized.

- Implant Detection
- TherapyGuide
- Capture Management (RA and RV)
- Sensing Assurance (RA and RV)
- Lead Monitor (RA and RV) with Auto Polarity Switch
- CareLink connectivity*

Minimizing Unnecessary RV Pacing

Promotes intrinsic conduction by reducing unnecessary RV pacing.

Search AV+ up to 600 ms

AT/AF Management

Pacing therapies to help manage atrial tachyarrhythmias and alleviate symptoms.

- ModeSwitch with Blanked Flutter Search
- Non-competitive Atrial Pacing (NCAP)

Additional Pacing Features

- Rate Hysteresis
- Sleep Function
- PVC Response
- Ventricular Safety Pacing (VSP)
- Dual Zone Rate Response Pacing with Rate Profile Optimization

Diagnostics

- Quick Look II
- Histogram Reports
- Atrial and Ventricular Episodes including EGMs
- Additional Clinician Selected Diagnostics



Model	SPDR01	SPDRL1 (Longevity)
M (g)	27.1	31.3
V (cc)	12.1	13.1
Size (mm) (HxWxD)	44.7x47.9x7.5	45.4x52.3x7.5
Connector	IS-1 BI/UNI	IS-1 BI/UNI

MRI SureScan

Full Body 1.5 and 3T MRI with any MRI SureScan lead:

- No MRI scan exclusion zone and no scan duration restriction
- MRI scan possible for the entire life of the system
- No Patient size restriction and no condition restrictions (e.g. fever)





^{*}Medtronic CareLink™ network are forward compatible; release date to be confirmed