# ATTESTA DR MRI TMI SURESCANTM

## PACEMAKERS (IPG)

## Completely Automatic - Simple to Use

Automatically monitors for implant detection and continuously adapts key device parameters to ensure the therapies are optimized.

- Implant Detection
- TherapyGuide
- Capture Management (RA and RV)
- Sensing Assurance (RA and RV)
- Lead Monitor (RA and RV) with Auto Polarity Switch
- CareLink connectivity\*

## **Minimizing Unnecessary RV Pacing**

Promotes intrinsic conduction by reducing unnecessary RV pacing.

- Managed Ventricular Pacing Mode (MVP): AAI(R) → DDD(R)
- Search AV+ up to 600ms

# AT/AF Management

Pacing therapies to help manage atrial tachyarrhythmias and alleviate symptoms.

- ModeSwitch with Blanked Flutter Search
- Post Mode Switch overdrive Pacing (PMOP)
- Atrial Preference Pacing (APP)
- Conducted AF Response (CAFR)
- Non-competitive Atrial Pacing (NCAP)

# **Additional Pacing Features**

- Rate Hysteresis
- Sleep Function
- PVC Response
- Ventricular Safety Pacing (VSP)
- Sinus Preference
- High Upper tracking rate up to 210 min<sup>-1</sup> for pediatric indications
- Rate Drop Response with 2 detection algorithms
- Dual Zone Rate Response Pacing with Rate Profile Optimization

## **Diagnostics**

- Quick Look II
- Cardiac Compass Trends
- Histogram Reports
- Atrial and Ventricular Episodes including EGMs
- Additional Clinician Selected Diagnostics



Model	ATDR01	ATDRS1 (Small)	ATDRL1 (Longevity)
M (g)	27.1	23.6	31.3
V (cc)	12.1	11.1	13.1
Size (mm) (HxWxD)	44.7 x 47.9 x 7.5	44.7 x 42.9 x 7.5	45.4 x 52.3 x 7.5
Connector	IS-1 BI/UNI	IS-1 BI/UNI	IS-1 BI/UNI

#### **MRI SureScan**

#### Full Body 1.5 and 3T MRI with any MRI SureScan lead:

- No MRI scan exclusion zone and no scan duration restriction
- MRI scan possible for the entire life of the system
- No Patient size restriction and no condition restrictions (e.g. fever)





<sup>\*</sup>Medtronic CareLink™ network are forward compatible; release date to be confirmed